Selectiva™ GUIDEWIRES

FILS-GUIDE Selectiva™

CABLES GUÍA Selectiva™

FILI-GUIDA Selectiva™

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Selectiva™-LEDETRÅDE

Selectiva™ VOERDRADEN

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VODIACE DRÔTY Selectiva™

NeoMetrics, Inc.
2605 Fernbrook Lane N, Suite J
Plymouth, MN 55447
Phone: +1 763.559.4440
Fax: +1 763.559.7676

European Representative:
Emergo Europe
Molenstraat 15
2513 BH, The Hague
The Netherlands
Phone: +31.70.345.8570
Fax: +31.70.346.7299

SP-3025 Rev. K

Selectiva™ GUIDEWIRES

Caution:
Federal law (USA) restricts this device for sale by or on the order of physician.

- Disposable, For Single Patient Use Only
- Ethylene Oxide (EO) Sterilization, Non-Pyrogenic
- Read all directions prior to use
- Store in a cool, dry place

Warning:

- This device is not intended for use in the coronary or cerebral vasculature.
- Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.
- Reprocessing may compromise the integrity of the device and / or lead to device failure.

Description:
Selectiva Guidewires are constructed of nickel titanium alloys, stainless steel, platinum/tungsten alloys, and feature a lubricious coating.

Indications:
To facilitate the placement of therapeutic devices for peripheral vascular interventional procedures. This device is not intended for use in the coronary or cerebral vasculature.

Directions:

1. Inspect the guidewire prior to use for tip shape, bends, kinks or coil separation. Do not use if damaged.
2. Using sterile technique, localize and puncture the vessel with a needle cannula.
3. Remove the needle, leaving the cannula in place.
4. Insert distal end of the guidewire through the cannula and into the vessel.
5. Remove cannula, leaving guidewire within the lumen of the vessel.
6. Pass the catheter over the guidewire within the lumen of the vessel under fluoroscope guidance to the desired position.
7. Carefully remove the guidewire from the catheter.

Cautions:

- Single use only, do not re-sterilize.
- Do not advance Selectiva Guidewires against resistance until the cause of the resistance has been determined.
- Excessive force against resistance may result in damage to guidewire and catheter or vessel perforation.
- Guidewires, by nature of their construction, will collect blood and other foreign matter in lumen. No type of cleaning will completely remove this material. Therefore, they are intended for single patient use only.
- NeoMetrics, Inc. does not recommend a particular technique for the use of this device. The steps contained in the preceding directions discuss the Seldinger Technique for percutaneous entry and are for information purposes only. Each physician should evaluate their appropriateness according to individual patient conditions and his or her medical training and experience.
- Do not use product if sterile package is damaged.

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